

## **AM Diagnostics**

### **All Test**

#### **Nasal**

#### **FAQ's**

##### ***What is the All Test?***

The All Test rapid antigen test is a single use test kit intended to detect the SARS-CoV-2 that causes COVID-19 by using a self-collected nasal swab specimen from symptomatic individuals.

##### ***Is the All Test approved for use in Australia?***

Yes, the All Test has been approved for use in Australia by the TGA. Details can be found on the TGA's website.

##### ***How long will it take to complete the test?***

Please allow time to complete your test without any distractions. Results are ready in 15 minutes. Please do not read results earlier than 15 minutes or more than 20 minutes.

##### ***Is the test easy to do?***

Yes, the test kit has been designed for use for the "Lay Person" meaning no medical knowledge required. The Gardian Self Check app will provide you with step-by-step instructions on how to complete your test.

##### ***Can children use the test?***

Yes, the test kit is suitable for children over the age of 2. Children must be tested by their parent or guardian. Dependents can be added to your account on the Gardian Self Check App. Please keep test kits out of children's reach.

##### ***My test is positive, what should I do?***

If you test positive, you must attend your nearest PCR centre and self-isolate. Please follow your local state/territory guidelines. Do not repeat the test hoping to receive a negative result.

##### ***My results are invalid, what should I do?***

An invalid result could be due to insufficient sampling or incorrect technique. The test kit could also be faulty, please repeat the test following the instructions provided.

##### ***I've accidentally spilt some of the buffer solution on my skin, what should I do?***

Please wash the area with plenty of water. If irritation occurs, please contact your healthcare provider.

##### ***The packaging of the test kit is already opened, is it still ok to use?***

No, please discard the test kit. Always remember to inspect the contents of the test kit prior to use.

##### ***How do I store the test kit?***

Please store test kits in a dry place at 2-30 degrees. Do not freeze and avoid excess moisture.

##### ***Can I reuse the test kit?***

No, all components of the test kits are designed for single use only.

***My test results are negative, but I still have symptoms, what should I do?***

If your test results are negative and you are still experiencing symptoms, you can carry out another test. If symptoms persist, please attend your local PCR testing centre, and follow your local state/territory guidelines.

***How accurate is the All Test kit?***

The All Test kits are 98% accurate, 95.9% sensitivity and 99.4% specificity.

***When is the best time to carry out the test?***

The best time to carry out the test is within the first 7 days of onset of symptoms.

***Can the All Test detect different strains of COVID-19?***

Yes, the All Test can detect the Alpha (UK strain), Delta (Indian), Gamma, VUI-21 ARP-03 (Indian strain) and Beta (South African strain).

***I have an issue with the test kit, who can I contact?***

If you have an issue with the performance or usability of the test kit, you can contact the TGA on 1800 809 361 or email [iris@tga.gov.au](mailto:iris@tga.gov.au)

***How do I dispose of the used test kit?***

Please dispose of used contents in the bio safety bag supplied.

***I've opened the test kit; do I need to use it straight away?***

Once the test cassette has been opened it should be used within the hour.